

PHOTOGRAPHY FOCUS




REFLECT

“These paediatric photography mirrors are the best length and shape on the current market for achieving the highest standard occlusal image. They fully reflect the light achieving sharp and detailed images.”

Zoë Kakiziba — Dental Photography Team Leader, Leeds Dental School


INTRA ORAL MIRRORS

With a solid chrome coating on 3 mm transparent glass, higher reflectivity is achieved – giving sharper, precise images. The mirrors withstand repeated steam and autoclavable sterilisation up to 137 °C. Choose from a range of six individually sold shapes.



SILICONE COATED MIRROR HANDLE

Adjustable mirror handle made from stainless steel, coated in Teflon silicone. Suitable for ultrasonic bath, washer and autoclave up to 137 °C.



DB05-0326T


Blue Handle
165 mm

DB05-0327

Pink Handle
165 mm

GLASS MIRRORS WITH HANDLE

These glass mirrors feature a stainless steel handle for ease of use. Autoclavable up to 137 °C.



DB04-0158

Occlusal
95 mm x 65 mm x 40 mm

DB04-0159

Buccal
115 mm x 40 mm

PHOTO MIRROR WITH HANDLE HAHNENKRATT

RHODIUM FRONT SURFACE + ULTRA FS SURFACE

The photo mirrors for intraoral recordings show an exactly sharp mirror image down to the smallest detail. Together with the brilliant natural colors and the pleasantly bright reflection, you get intraoral images of great significance.

The photo mirrors are available in different models, sizes and as FS Rhodium and FS ULTRA versions.




CONTRAST

DB contrastors improve your anterior, occlusal and lateral photographs by eliminating unwanted details such as the tongue, nostrils, lips and moustaches. They feature an exclusive anatomical format and opaque finish; and are autoclavable up to 137 °C.

CONTRASTORS

These one piece contrastors allow for anterior and occlusal photographs and can also be used for lateral photographs. Autoclavable up to 137 °C



DB04-0157

Palatal Contrastor
240 mm x 65 mm



Palatal Contrastor
240 mm x 55 mm

RETRACT

“I've tried out the DB04-0154A stainless steel retractor a few times now and can give nothing but really positive feedback from myself, my staff and my patients. It is easy to use, brilliant for the autoclave, and the slim-line feel allows for maximum patient comfort when positioning it. In addition, the handle design reduces any risk of your hand slipping when retracting, I give this product an all-round thumbs up!”

Ailbhe McMullin — Consultant Orthodontist, Royal Manchester Children's Hospital

PHOTO LIP RETRACTOR STAINLESS STEEL

The durability of DB stainless steel cheek retractors make them a high quality, yet cost effective solution. Featuring a satin finish to prevent glare, these retractors are autoclavable up to 137 °C.

The slim line design allows maximum exposure for occlusal photographs, and easy insert for patient comfort.




DB04-0154 A

Photo Lip Retractor

PHOTO CHEEK RETRACTOR STAINLESS STEEL

One large and one small end. Provides excellent visibility of the entire dental arch.

Fully autoclavable at 134–137 °C.



DB04-0199

Stainless Steel
Cheek Retractor
110 mm x 50 mm x 40 mm

DB04-0192

Stainless Steel
Cheek Retractor
170 mm x 40 mm x 35 mm

PHOTO LIP RETRACTORS

Innovative design provides maximum exposure to occlusal photographs. Autoclavable up to 121 °C. Sold as each.



DB04-0154
Photo Lip Retractor



05047
K&K Occlusal
Retractor Large




05046
K&K Occlusal
Retractor Small

LONG HANDLED CHEEK RETRACTOR

Autoclavable up to 121 °C Sold in pairs.



DB04-0190
Long Handled
Cheek Retractor
110 mm x 35 mm



DB04-0191
Long Handled
Cheek Retractor
115 mm x 45 mm

GOLD SERIES CHEEK RETRACTORS

Introducing our new Gold Series cheek retractors which have been designed and manufactured in the UK utilising a revolutionary new medical grade material which will withstand multiple high heat sterilisation cycles up to 137 °C without deformation.

Initially available in two of our best-selling designs these new products are guaranteed to impress! Designed to give a better buccal view. Sold in pairs.



DB04-0174 A

Small Gold Series
Cheek Retractors
145 mm x 40 mm x 35 mm

DB04-0175 A

Large Gold Series
Cheek Retractors
175 mm x 40 mm x 35 mm

DB04-0178

Small Gold Series
Double Ended V Retractors
160 mm x 35 mm x 35 mm

DB04-0178L

Large Gold Series
Double Ended V Retractors
175 mm x 55 mm x 35 mm

BEST SELLING RETRACTORS

SOFT SPANDEX CHEEK EXPANDERS

The soft material of the soft spandex cheek expander features a matt finish to reduce glare, and allows for greater flexibility; reducing the fracture risk associated with harder plastic models. Autoclavable up to 137 °C.



DB04-0186
Small



DB04-0185
Large

SMALL CHEEK RETRACTOR


Sized for adult and child. Durable plastic design. Autoclavable up to 137 °C.



DB04-0183
100 mm x 45 mm x 40 mm


PHOTO CASSETTES

The soft material of the soft spandex cheek expander features a matt finish to reduce glare, and allows for greater flexibility; reducing the fracture risk associated with harder plastic models. Autoclavable up to 137 °C.



DB06-5028

Medium low profile cassettes,
holds 3 mirrors securely
198 mm x 165 mm x 37 mm




DB06-5029

Large low profile cassettes,
holds 4 mirrors and 2 cheek retractors
290 mm x 183 mm x 37 mm


INDIVIDUAL PHOTO CASSETTES

DB Cassettes protect your photographic mirrors in the autoclave. Made of stainless steel with silicone slots, the mirrors are kept securely in place with no movement – reducing the risk of damage.

Available in 3 sizes.



DB06-5029 S
140 mm x 80 mm x 15 mm



DB06-5029 M
170 mm x 70 mm x 15 mm



DB06-5029 L
200 mm x 80 mm x 15 mm



MEDEQ

MEDICAL EQUIPMENT AND SUPPLIES

INFORMATION AND ORDERS

info@medeq.ee
www.medeq.ee

+372 622 8111
+372 5693 5646