

GC



MEDEQ
MEDICAL EQUIPMENT AND SUPPLIES

GC Orthodontics
breketid, tuubid ja kaared

Iseligeeruvad breketid

EXPERIENCE

SOODNE

ESTEETILINE

ESTEETILINE

EXPERIENCE metal

SL metallbreket

- interactive SL- system
- low-friction/control combination



- MBT, Roth, Extra Torque, Andrews, GET
- .018", .022"
- medium ja mini suuruses

EXPERIENCE metal Rhodium Coated

Roodiumiga kaetud SL breket

- interactive SL- system
- low-friction/control combination
- Rhodium coated
- Rhodium clip



- MBT, Roth, GET
- .018", .022"
- medium ja mini suuruses

EXPERIENCE ceramic

Keraamiline SL breket

- interactive SL- system
- low-friction/control combination



- MBT, Roth, Extra Torque, GET
- medium ja mini suuruses

Tuubid


EXPERIENCE



Experience SL tuubid

- iseligeeruvad 6ndate tuubid
- MBT, Roth, Extra Torque
- .018", .022"
- single ja double tuubid

LEGEND



Legend tuubid

- konventsionaalsed 6-ndate ja 7-ndate tuubid
- LP- low profile ja LE-large entrance variandid
- MBT, Roth, Andrews, GET
- .018", .022"



GC Ortho Connect™

Metall- ja keraamiliste breketite ning tuubide kleepimiseks

- sisaldab primerit
- valguskõvastuv

Konventsionaalsed breketid

SOODNE

PALJUDE PRESCRIPTIONITEGA

ESTEETILINE

AXCESS


Metallbreket
AXCESS



- MBT, Roth
- .018", .022"

LEGEND


Metallbreket
LEGEND



- MBT, Roth, Alexander, Extra Torque, Bio- Progressive, GET
- .018", .022"
- medium ja mini suuruses

LEGEND Rhodium Coated

**Roodiumiga kaetud
metallbreke**
LEGEND RC



- Rhodium coated
- MBT, Roth
- 018", .022"
- medium ja mini suuruses

ESTEETILINE

Chíc.


**Keraamiline breket
CHIC**



- MBT, Roth
- .018", .022"
- medium ja mini suurus

Tuubid

LEGEND



Legend tuubid

Konventsionaalsed 6-ndate
ja 7-ndate tuubid

- LP- low profile ja LE-large
entrance variandid
- MBT, Roth, Andrews, GET
- .018", .022"



GC Ortho Connect™

Metall- ja keraamiliste breketite
ning tuubide kleepimiseks

- sisaldab primerit
- valguskõvastuv

Kaared

Ümmargune thermokaar INITIALLOY/ INITIALLOY RHODIUM

- saadaval Rhodium coated
- lightest forces
- superelastic
- .012"– .018"



Kandiline thermokaar BIO-EDGE

- shape memory
- .016" x .016" – .021" x .028"



TMA (Beta Titanium)

- .016" – .021" x .028"



Kandiline thermokaar BIO-ACTIVE BIO-ACTIVE LIGHT BIO-ACTIVE RHODIUM

- saadaval Rhodium coated
- shape memory
- progressive forces
- .016" x .016" – .021" x .028"



8-stranded SS MULTIFLEX

- .016" x .022"– .021" x .025"



Monostranded SS/ Monostranded SS RHODIUM

- saadaval Rhodium coated
- .014"– .021"– .025"



Available Prescriptions by Brackets Type

Brackets/Prescription	GET	Ext. Torque	Roth Type	MBT Type	ACO	And. Type	Bio. Prog.	Bi Dim.
Experience mini metal								
Experience mini metal Rh								
Experience metal								
Experience metal Rh								
Experience mini ceramic								
Experience ceramic								
Chic.								
Chic. mini								
LEGEND mini*								
LEGEND mini Rh								
LEGEND medium**								
LEGEND medium Rh								
Axess								
Legacy								
Legacy mini								

*LEGEND mini is also available in Planche Type, V-Slot Roth Type, Aeksander Type


**LEGEND medium is also available in Bio-Progressive V-Slot, Bio_Dimentional V-Slo, Polar

Angulation and Rotation Information by Prescription

Maxillary	U1			U2			U3			U4			U5			U6			U7		
Prescription	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.
GET	12°	5°	-	8°	9°	-	0°	11°	-	-10°	4°	2° _d	-10°	4°	2° _d	-20°	0°	10° _d	-25°	0°	6° _d
Extra Torque	18°	5°	-	10°	8°	-	0°	8°	2° _m	-10°	4°	2° _d	-10°	4°	2° _d	-20°	0°	20° _d	-20°	0°	20° _d
Roth Type	12°	5°	-	8°	9°	-	-2°	13°	-	-7°	0°	-	-7°	0°	-	-14°	0°	14° _d	-14°	0°	14° _d
MBT Type	17°	5°	-	10°	9°	-	-7° 0	8°	-	-7°	0°	-	-7°	0°	-	-14°	0°	10° _d	-19°	0°	8° _d
ACO	12°	5°	-	10°	9°	-	-7°	10°	2° _m	-9°	0°	-	-9°	0°	-	-14° -20°	0°	-10° -14°	-14° -20°	-	10° _d 14° _d 10° _d
Andrews Type	7°	5°	-	3°	9°	-	-7°	11°	-	-7°	2°	-	-7°	2°	-	-9°	5°	10° _d	-9°	5°	10° _d
Bioprogressive Type	22°	5°	-	14°	8°	-	7°	10°		-7°	0°		-7°	0°		-10°	0°	15° _d	-10°	0°	12° _d

Angulation and Rotation Information by Prescription

Mandibular	U1			U2			U3			U4			U5			U6			U7		
Prescription	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.	T.	A.	R.
GET	-1°	2°	-	-1°	2°	-	-3°	3°	-	-17°	1°	4° _d	-22°	1°	4° _d	-30°	0°	4° _d	-30°	0°	4° _d
Extra Torque	-1°	0°	-	1°	0°	-	-3°	7°	2° _m	-7°	-1°	2° _d	-7°	-1°	2° _d	-10° _d -14° _d	0°	4° _d	-10°	0°	-
Roth Type	-1°	2°	-	-1°	2°	-	-11°	7°	-	-17°	0°	-	-22°	-1°	-	-30°	-1°	4° _d	-30°	-1°	4° _d
MBT Type	-6°	0°	-	-6°	0°	-	-6° 0	3°	-	-12°	2°	-	-17°	2°	-	-20°	0°	-	-10°	0°	-
ACO	-6°	0°	-	-6°	0°	-	-8°	3°	2° _m	-12°	2°	-	-17°	-1°	4° _d	-25° _d -20°	-1°	-	-20° -25°	-1°	-
Andrews Type	-6°	2°	-	-6°	2°	-	-11°	5°	-	-17°	2°	-	-22°	2°	-	-30°	2°	-	-31°	2°	-
Bioprogressive Type	-1°	0°	-	-1°	0°	-	7°	5°		-11°	0°	-	-17°	0°		-27°	0°	5° _d	-27°	0°	6° _d



'GC'

GC Orthodontics ametlik edasimüüja
Eestis on Medeq OÜ

Rohkem infot ja tellimine:

MEDEQ
MEDICAL EQUIPMENT AND SUPPLIES

Medeq OÜ
+372 622 8111
info@medeq.ee
www.medeq.ee