

Oral Screens

STOPPI® (for children from the age of 2 to 4)



Stoppi® is a dental-care alternative to existing sucking habits. STOPPI® has no foreign bodies between the incisors. Therefore the upper and lower jaw can develop normally and the teeth can find their correct position without being disturbed.

DESCRIPTIO	ORDER NUMBER
N Unit	EC M0166
Kit - 6 pcs	EC M0165

Oral Screen – Standard

Even if teeth do not suffer from major malpositioning damaging thumb-sucking habits should be corrected with the assistance of the Oral screen.

The Oral screen:

- promotes the self-healing of existing anomalies
- helps eliminate speech defects
- corrects malfunctions of the tongue
- eliminates harmful mouth breathing
- leads to normal lip closure
- strengthens lip muscles
- supports myofunctional interventions

DESCRIPTION	ORDER NUMBER
Small model	EC M0150
Large model	EC M0151
Soft small model (clear)	EC M0160
Soft large model (clear)	ECM0161

Oral Screen with Bite Caps

Thumb-sucking frequently produces large overbite of the anterior teeth, or deteriorates an existing retroposition of the mandible. In this case the Oral screen with a bite cap is indicated. The bite cap prevents the Oral screen from tipping. The mandible is simultaneously guided forward.

DESCRIPTION	ORDER NUMBER
Small model	EC M0152
Large model	EC M0153
Soft small model (clear)	EC M0162
Soft large model (clear)	EC M0163

Oral Screen with Wire Guard

This form of Oral screen is indicated for an anterior open bite caused by the thumb or by frequent use of a pacifier. The tongue crib prevents the tongue from inserting between the dental arches and helps closing the open bite.

DESCRIPTIO	ORDER NUMBER
N Small model	EC M0154
Large model	EC M0155

Oral Screen with Pearl

The pearl, held by a wire frame. Guides a tongue with weak muscles toward the palate. This special Oral screen is indicated in the case of Downs syndrome, speech defects or cerebral palsy.

DESCRIPTIO	ORDER NUMBER
N Small model	EC M0156
Large model	ECM01537



Ordering and information:

Medeq OÜ
info@medeq.ee

+372 622 8111
www.medeq.ee