

PHOTOGRAPHY FOCUS & ALL YOU NEED FOR ALIGNERS

price list 2025 & 2026 without VAT

ORTHOQ





190€



190€



190€

ixionTM
Instruments

IXION ALIGNER PLIER **TORQUE**

IX887

Creates a horizontal indentation, 0.75mm x 1.5mm, for individual root torque.

IXION ALIGNER PLIER **ROTATION**

IX888

Creates a vertical indentation, 1.5mm x 0.75mm, for individual tooth rotation.

IXION ALIGNER PLIER **RETENTION**

IX889

Creates a circular indentation, 1mm in diameter, for increased retention in the retainer and aligner system.



90€

IXION ALIGNER PLIER SCISSORS

IX892

These lightweight scissors are used for removing excess thermoformed material, and for the smooth trimming and shaping of aligners. Their streamline design, graduated curve and tapered beak makes for easy use - even in tight areas.



220€

IXION ALIGNER PLIER TEAR DROP

IX890

Creates a hook in the gingival margin for the use of elastomerics.



220€

IXION ALIGNER PLIER HOLE PUNCH

IX891

Creates a circular cutout to allow bondable attachments to be used.

IXION ALIGNER PLIER INSTRUMENT SET

IX886

The popular choice for premium instrumentation of orthodontists globally. This kit includes all 5 of our Ixion Aligner Pliers in a metal Ixion branded presentation tin.

959€



www.orthoq.ee

DB ORTHODONTICS **INTERPROXIMAL STRIPPING (IPR)**

An innovative Interproximal Reduction system that makes interproximal enamel removal / stripping safer and more accurate without creating sharp corners or subgingival ledges.



Airo Strips Intro pack
DB05-0908

145€

Designed for use with all orthodontic treatments, these strips curve and fly along the natural contours of the teeth to avoid cutting lips and gums, prevent excessive enamel removal and to maximize patient comfort.



Yellow 0.07mm Single-Sided

DB05-0900

Red 0.12mm Double-Sided

DB05-0901

Dark Blue 0.15mm Double-Sided

DB05-0902

Green 0.20mm Double-Sided

DB05-0903

Clear 0.10mm Single-Sided

DB05-0904

Cyan 0.15mm Single-Sided

DB05-0905

Purple 0.25mm Double-Sided

DB05-0906

Brown 0.30mm Double-Sided

DB05-0907



45€

INCREMENTAL **THICKNESS GAUGE**

DB05-0915

Made from Stainless Steel - 0.1, 0.2, 0.25, 0.3, 0.4, 0.5 mm. Measure the exact amount of interproximal reduction achieved.

DENTASONIC BASIC SET ALIGNER

DS64A-204-2Y

Basic introduction IPR set with strips, gauge.
Water cooled Hand Piece.

Content:

- DSS-01, 1x DSM-02, 1x DSM-03, 1x DSM-04, 1x RDS-05, 2 sides, 3X RPS-15
- 1x DSG-204 Twist water-cooling HP, 4:1, STROKE 1,1mm 1x MG-06,
- Measuring gauge with 6 thicknesses: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5mm
- 1 extraction tip + 1 Spray nozzle
- 1x Water hole needle, 1x HP Instruction booklet, 1x Worldwide Warranty 1x IFU, Instructions for use



750€

DENTASONIC PROFESSIONAL KIT 2

DSPK2-2Y

Professional Kit 2 with Twist water-cooling HP, 1 additional reciprocating head, 2 years warranty and reparation voucher.

The kit includes:

- 1x DSS-01, DSS-02, DSM-0.3, DSM-0.4, RDS-04, RDS-0.5mm 2 sides
- 1x DSS-005, DSM-01, DSM-02, DSM-03, RDS-009L, RDS-010Lmm 1 side 3x RPS-9, 3X RPS-15
- 1x DSG-204 Twist water-cooling HP, 4:1, STROKE 1,1mm
- 1x DSH-10 reciprocating Head
- 2x DB-R20 Cup, 2x DB-R21 Point, Debonding polishers
- 1x MG-06, Measuring gauge with 6 thicknesses: 0.1, 0.15, 0.2, 0.3, 0.4, 0.5mm
- 1x ST-5, Stainless steel SteryStrip Holder
- 1x MG-1 extraction tip + 1x SN-1 Spray nozzle, 1x WN-1 Water hole needle
- 2x CS-60-40, 1x CS-20-7, Finishing CombiStrip
- 1x HP Instruction booklet, 1x Worldwide Warranty, 1x Leaflet 1x IFU, Instructions for use.



995€

CERAMIC BUTTON FREDDY™ CURVED BASE

44-1010-0000

The Freddy™ ceramic button is the perfect clear utility button for all aesthetic appliances.

10 pcs per package.

- Works lingually as well as labially
- Mesh structure of base makes a strong bonding
- Bonds well with most common adhesives
- Anatomy and size base fits all teeth
- Enough space around tie wing for elastics
- Maintains the transparent look with aligners as anchorage for elastics



105€

LINGUAL BUTTON AND LINGUAL LOOP

Bondable lingual buttons are used in dental movements such as retractions, rotations, pull and anchoring systems. This retains the elastic, wires, chains and ligatures. 10 pcs per package.



27€

Lingual button curved base

42-0010-0000

Lingual button flat base

42-0020-0000

Lingual loop for curved teeth

42-0030-0000

Lingual loop for flat teeth

42-0040-0000





97€

GC ALIGNER CONNECT

90-4100-0000

GC Aligner Connect was specifically developed for orthodontic aligner attachments. GC Aligner Connect features enhanced aesthetics that imitate enamel color shades and high durability, which limits long term wear or potential detachment during treatment.



GC ORTHO CONNECT

90-1100-0000

A new orthodontic light-cured adhesive for brackets and bondable buccal tubes.

GC ORTHO CONNECT FLOW

90-3100-0000

Light cure orthodontic adhesive for lingual retainers, stops and aligner attachments - easy to handle. GC Ortho Connect Flow's **primer is integrated** into the paste, allowing to save valuable time in the bonding process.

ORCHESTRA INTRA AND EXTRAORAL **ELASTICS**

Our range of intra and extraoral elastics are made from the highest quality medical grade latex. The stated force corresponds to the force three times its original listed diameter.

Intraoral elastics are packaged by 5000 pcs.

Extraoral elastics are packaged by 2500 pcs.

Each bag contains 1 insertion plastic instrument.

Latex free elastics are listed on the chart below.

45€



INTRAORAL				EXTRAORAL
SIZE/FORCE	2.5 oz/ 70 g	4.0 oz/ 113 g	6.0 oz/170 g	8.5 oz/ 240 g
1/8" / 3mm	Harp 85-0321-0000	Xylophone 85-0331-0000	Tuba 85-0341-0000	-
3/16" / 4mm	Tambourine 85-0421-0000	Piano 85-0431-0000	Trumpet 85-0441-0000	Trumpet 85-0441-0000
3/16" / 4mm Latex Free	-	Jazz 85-0430-0000	Rock-n-Roll 85-0440-0000	-
1/4" / 6mm	Oboe 85-0621-0000	Flute 85-0631-0000	Baritone 85-0641-0000	French Horn 85-0661-0000
1/4" / 6mm Latex Free	-	Hip-Hop 85-0630-0000	Classical 85-0640-0000	-
5/16" / 8mm	Cello 85-0821-0000	Clarinet 85-0831-0000	Saxophone 85-0841-0000	Guitar 85-0861-0000
3/8" / 10mm	Cymbals 85-1021-0000	Tympani 85-1031-0000	-	Violin 85-1061-0000

GC ELASTOMERIC CHAINS

Colour	Short	Medium	Long
Clear	86 - 1010 - 0000	86 - 1020 - 0000	86 - 1030 - 0000
Gray	86 - 2010 - 0000	86 - 2020 - 0000	86 - 2030 - 0000

Made from pure high-density and high-purity materials, this chain is precision-cut and the intervals between links are strictly identical. GC Ortho Chain exerts light and continuous forces for optimal tooth movement.

Available in 2m spools grey or clear.



Short

Medium

Large

29€



MULTI-FA

Multi functional Appliance for **CL II and CL III Correction**

Simple, Precise, Adaptable, Efficient

One appliance for **90% cases**

175€

Available size:

Multi- Fa **Regular** (18-28 mm)


92-5001-0000

Multi- Fa **Long** (23-36 mm)

92-5003-0000

Short (12-17 mm)

92-5002-0000



Watch Dr. Mazor's videos
about Multi-Fa here.



ORTHO EASE ALIGNER CHEWIES

DB04-4231 Plain

DB04-4230 Mint

Chewies help the patient seat the aligner into place securely and comfortably. They also stimulate the periodontal ligaments to aid effective tooth movement. Each package contains 10 boxes of 2 chewies in a hygienic sealed bag.

ORTHO EASE ALIGNER RETRIEVERS

DB04-4220 White

Remove aligners and retainers quickly and comfortably with easy to use retrievers. Each package contains 10 retrievers in a hygienic sealed bag



3.5€

5.5€

27€



27€



MYOCLEAN ANTIBACTERIAL CLEANING TABLETS FOR DENTAL APPLIANCES

301212- Package of 14 tablets

301224- Package of 24 tablets

Mouthguards are a perfect breeding ground for bacterial growth. Bacteria get trapped in crevices and pores of the material and if not cleaned properly, the germs can flourish and multiply.

Myoclean™ keeps your mouthguard hygienically clean and removes stains, while leaving it with a fresh minty taste.

INSTRUCTIONS FOR USE: Add 1 tablet to a glass of lukewarm water. Leave the appliance to soak for a minimum of 15 minutes.

AUTOGRAPH WHITENING

Autograph Whitening is a professional teeth whitening system provided to you by your dentist or orthodontist.

Safe, gentle, effective, and easy to use – our whitening system delivers results within the first two weeks. The product contains **16% carbamide peroxide** and requires only **1 hour of wear time** per session.



35€

AG02-0004

2 x 3ml Syringe Top-Up Pack
(16% Carbamide Peroxide)

65€


AG02-0003

4 x 3ml Syringe Patient Kits
(16% Carbamide Peroxide)

INTRA ORAL MIRRORS

89€

With a solid chrome coating on 3 mm transparent glass, higher reflectivity is achieved - giving sharper, precise images. The mirrors withstand repeated steam and autoclavable sterilisation up to 137 °C. Choose from a range of six individually sold shapes.



DB04-0160
Palatal Adult
70 mm x 100 mm

DB04-0161
Palatal Child
60 mm x 85 mm

DB04-0161/C
Palatal Child/Adult
56 mm x 115 mm x 65 mm

DB04-0162
Buccal
40 mm x 150 mm

DB04-0163
Lingual
60 mm x 150 mm

DB04-0164
Large Palatal
70 mm x 135 mm

PHOTO MIRROR WITH HANDLE HAHNENKRATT

RHODIUM FRONT SURFACE

The photo mirrors for intraoral recordings show an exactly sharp mirror image down to the smallest detail. Together with the brilliant natural colors and the pleasantly bright reflection, you get intraoral images of great significance.



Rhodium surface with bent handle

190€

CONTRASTORS

89€

These one piece contrastors allow for anterior and occlusal photographs and can also be used for lateral photographs. Autoclavable up to 137 °C. Sold in pairs.

240 mm



DB04-0157

Palatal Contrastor
240 mm x 65 mm

240 mm



Palatal Contrastor
240 mm x 55 mm

PHOTO LIP RETRACTORS

Innovative design provides maximum exposure to occlusal photographs. Autoclavable up to 121 °C. Sold as each.



39€

RET-OCCL-SMALL

K&K Occlusal
Retractor Small

RET-OCCL-MED

K&K Occlusal
Retractor Large



35€

DB04-0154

Photo Lip Retractor

SMALL CHEEK RETRACTOR

Sized for adult and child. Durable plastic design. Autoclavable up to 137 °C.



DB04-0183

100 mm x 45 mm x 40 mm

29€

GOLD SERIES CHEEK RETRACTORS

Introducing our new Gold Series cheek retractors which have been designed and manufactured in the UK utilising a revolutionary new medical grade material which will withstand multiple high heat sterilisation cycles up to 137 °C without deformation.

Initially available in two of our best-selling designs these new products are guaranteed to impress! Designed to give a better buccal view. Sold in pairs.



DB04-0174 A
Small Gold Series
Cheek Retractors
145 mm x 40 mm x 35 mm

DB04-0175 A
Large Gold Series
Cheek Retractors
175 mm x 40 mm x 35 mm

DB04-0178
Small Gold Series
Double Ended V Retractors
160 mm x 35 mm x 35 mm

DB04-0178L
Large Gold Series
Double Ended V Retractors
175 mm x 55 mm x 35 mm

**BEST SELLING
RETRACTORS**

PHOTO CASSETTES


The soft material of the soft spandex cheek expander features a matt finish to reduce glare, and allows for greater flexibility; reducing the fracture risk associated with harder plastic models. Autoclavable up to 137 °C.



DB06-5028

Medium low profile cassettes,
holds 3 mirrors securely
198 mm x 165 mm x 37 mm

180€



DB06-5029

Large low profile cassettes,
holds 4 mirrors and 2 cheek retractors
290 mm x 183 mm x 37 mm

235€



Follow and tag us on Instagram
@orthoq_baltics



INFORMATION AND ORDERS

info@orthoq.ee
+372 5693 5646
+372 5662 6750
www.orthoq.ee

