

Watch Dr. Mazor's videos
about Multi-Fa here.



ALL YOU NEED FOR MULTI-FA TREATMENT

Multi functional for CL II and CL III Correction

2026 price list without VAT



ORTHOQ

MULTI-FA PROTOCOL



MULTI-FA

Multi functional Appliance for **CL II and CL III Correction**
Simple, Precise, Adaptable, Efficient
One appliance for **90% cases**

175€

Available size:

Multi- Fa **Regular** (18-28 mm)


92-5001-0000

Multi- Fa **Long** (23-36 mm)

92-5003-0000

Short (12-17 mm)

92-5002-0000




Multi-Fa stops: **MULTI-FA-STOPS** (big stops for spring activation) **20€**

Screwdriver: **MULTI-FA-STOP-HANDLE** (Multi-Fa handle for stops) **10€**

Outer spring: **MULTI-FA-SPRING** (Multi-Fa NiTi coil spring 0.18 or 0.16) **20€**

Inner spring for opening gaps: **61-0140-0013** (Open spring 91cm / 3' spool; .010"x.030") **15€**

Multi Fa for molar tubes: **92-5004-0000** (Multi Fa appliance for molar tubes or molar bands) **175€**



CERAMIC BUTTON FREDDY™ CURVED BASE

44-1010-0000

The Freddy™ ceramic button is the perfect clear utility button for all aesthetic appliances.

10 pcs per package.

- Works lingually as well as labially
- Mesh structure of base makes a strong bonding
- Bonds well with most common adhesives
- Anatomy and size base fits all teeth
- Enough space around tie wing for elastics
- Maintains the transparent look with aligners as anchorage for elastics



105€

LINGUAL BUTTON AND LINGUAL LOOP

Bondable lingual buttons are used in dental movements such as retractions, rotations, pull and anchoring systems. This retains the elastic, wires, chains and ligatures. 10 pcs per package.

27€



Lingual button curved base

42-0010-0000

Lingual button flat base


42-0020-0000

Lingual loop for curved teeth

42-0030-0000

Lingual loop for flat teeth

42-0040-0000



97€



GC ORTHO CONNECT

90-1100-0000

A new orthodontic light-cured adhesive for brackets and bondable buccal tubes.

For attaching Multi-FA to the tooth surface.




GC ORTHO CONNECT FLOW

90-3100-0000

Light cure orthodontic adhesive for lingual retainers, stops and aligner attachments - easy to handle. GC Ortho Connect Flow's **primer is integrated** into the paste, allowing to save valuable time in the bonding process.

For locking the Multi-FA telescope and joint.



ORCHESTRA INTRA AND EXTRAORAL **ELASTICS**

Our range of intra and extraoral elastics are made from the highest quality medical grade latex. The stated force corresponds to the force three times its original listed diameter.

Intraoral elastics are packaged by 5000 pcs.

Extraoral elastics are packaged by 2500 pcs.

Each bag contains 1 insertion plastic instrument.

Latex free elastics are listed on the chart below.

45€



INTRAORAL				EXTRAORAL
SIZE/FORCE	2.5 oz/ 70 g	4.0 oz/ 113 g	6.0 oz/170 g	8.5 oz/ 240 g
1/8" / 3mm	Harp 85-0321-0000	Xylophone 85-0331-0000	Tuba 85-0341-0000	-
3/16" / 4mm	Tambourine 85-0421-0000	Piano 85-0431-0000	Trumpet 85-0441-0000	Trombone 85-0461-0000
3/16" / 4mm Latex Free	-	Jazz 85-0430-0000	Rock-n-Roll 85-0440-0000	-
1/4" / 6mm	Oboe 85-0621-0000	Flute 85-0631-0000	Baritone 85-0641-0000	French Horn 85-0661-0000
1/4" / 6mm Latex Free	-	Hip-Hop 85-0630-0000	Classical 85-0640-0000	-
5/16" / 8mm	Cello 85-0821-0000	Clarinet 85-0831-0000	Saxophone 85-0841-0000	Guitar 85-0861-0000
3/8" / 10mm	Cymbals 85-1021-0000	Tympani 85-1031-0000	-	Violin 85-1061-0000

GC ELASTOMERIC CHAINS

Colour	Short	Medium	Long
Clear	86 - 1010 - 0000	86 - 1020 - 0000	86 - 1030 - 0000
Gray	86 - 2010 - 0000	86 - 2020 - 0000	86 - 2030 - 0000

Made from pure high-density and high-purity materials, this chain is precision-cut and the intervals between links are strictly identical. GC Ortho Chain exerts light and continuous forces for optimal tooth movement.

Available in 2m spools grey or clear.



Short

Medium

Large

29€





Follow and tag us on Instagram
@orthoq_baltics



INFORMATION AND ORDERS

info@orthoq.ee
+372 5693 5646
+372 5662 6750
www.orthoq.ee



Dr. Mazor's videos about Multi-Fa



Multi- FA protocol